

For higher density arrays, (approximately 13 μ spot size), fluorescent signal from 3 μ resolution scanning typically provides 16 pixels per feature which translates to improved confidence in the resulting data and increased number of reference CNVs detected when compared to 5 μ resolution scanning. When scanning the same arrays at 2 μ resolution, at least 36 pixels per feature are generated; which greatly enhances the fluorescent data leading to more accurate calls and increased confidence in the final examination and reporting of chromosomal aberrations.

Ordering information



For more information on testing services or to order stem cells from WiCell®, go to: <http://www.WiCell.org/>



For more information on high resolution scanning or to order an INNOPSYS InnoScan® microarray scanner go to: <http://www.innopsys.com/>

Authors and Contributors

Karen D. Montgomery, Ph.D, WiCell Research Institute, kmontgomery@wicell.org
 Seth Taapken, MS-Biotech, WiCell® Research Institute, staapken@wicell.org
 Ben Nisler, BS, WiCell® Research Institute, bnisler@wicell.org
 Thomas E. Strader, MS-Biotech, Heartland Biotech, tom.strader@heartlandbiotech.com

References


1. Stem cell lines which were formerly available for federally funded research need to be re-certified to ensure compliance with current regulations. Cell line WAO1 (H1) is an approved cell line available from WiCell®. See http://grants.nih.gov/stem_cells/registry/current.htm?id=20 for details. For more information, refer to: <http://stemcells.nih.gov/>
2. Genome Technology Magazine, Genetic Variation Technical Guide, February 2010

Disclaimers

WiCell® is a registered trademark of the WiCell Research Institute.
 InnoScan® is a registered trademark of INNOPSYS.
 Genepix® is a registered trademark of Molecular Devices
 CCHfusion® is a registered trademark of InfoQuant.
 NimbleGen® is a registered trademark of Roche.
 DNeasy® is a registered trademark of QIAGEN.
 SpeedVac® and NanoDrop® are registered trademarks of Thermo Scientific.

Applications note

Characterization of Stem Cell Copy Number Variation by aCGH at the Wisconsin International Stem Cell Bank



Fast Accurate Results Accelerate Stem Cell Science

Feature Extraction and Data Analysis

Feature extraction and segmentation were performed with NimbleScan® v2.5 software (Roche NimbleGen). Data analysis was performed with CGH Fusion® software (InfoQuant).

Typical aberration detection settings:

Algorithm: RBSv1.0 (robust circular binary segmentation)

p-value: 0.0010 or 0.010

Min. aberration length (probes): 5

Min. aberration length (Mbp): 0.0

Average log-ratio threshold: 0.0 or 0.2

Ratio plots were generated using CGH Fusion® software. A ratio plot is a visual representation of differences in the presence of hybridization to the microarray of the test sample and the reference sample. If there are equal amounts of test and reference sample, the dot is blue. If there is more test than reference (amplification), the dot is green. If there is more reference than test (deletion), the dot is red.




Figure 5. Comparison of ratio plots of the same microarray scanned at 5 μ (top), 3 μ (middle) and 2 μ (bottom) resolution on the InnoScan®900AL. Note, for example, that in chromosome 6 (in box) successively more aberrations (in green or red) were detected at 3 μ resolution than at 5 μ resolution; and still more were detected at 2 μ resolution than at 3 μ resolution.




Table 1. Flow-chart of protocol for aCGH sample preparation and analysis

Procedure

Cells are removed from cryogenic storage, allowed to come to room temperature, and propagated in a 6-well plate for 5-10 days per WiCell® Standard Operating Procedure SOP-CC-005. DNA is extracted using DNeasy® Blood & Tissue Kit (QIAGEN) per WiCell® SOP-CH-014. DNA quality and quantity are measured using the Nanodrop-1000 spectrophotometer (Thermo Scientific). Sample integrity and absence of RNA contamination is measured using DNA gel electrophoresis and imaging per SOP-CH-007. Samples are then labeled with Cy-3 and Cy-5 fluorescent dyes (Tri-Link) followed by overnight Klenow (New England Biolabs) amplification. DNA is isopropanol precipitated, dried in a SpeedVac® (Thermo Scientific) and reconstituted. Alexa Fluor® (Life Technologies) dyes are added followed by another isopropanol precipitation and SpeedVac® centrifugation until dry. Dried pellets are then rehydrated in VWR water.

Hybridization solution, including alignment oligo, is made, sample quantified using the nucleic acid feature on the NanoDrop® (Thermo-Fisher) spectrophotometer and the mixture adjusted to ensure optimal concentration for hybridization. Equal amounts of test and reference sample are combined in VWR water; after which the appropriate amount of hybridization solution master mix is added for a total volume of 18 μ L (385K array) or 20 μ L (3x720K array).

Samples are then labeled with Cy-3 and Cy-5 fluorescent dyes (Tri-Link) followed by overnight Klenow (New England Biolabs) amplification. DNA is isopropanol precipitated, dried in a SpeedVac® (Thermo Scientific) and reconstituted. Alexa Fluor® (Life Technologies) dyes are added followed by another isopropanol precipitation and SpeedVac® centrifugation until dry. Dried pellets are then rehydrated in VWR water.

Hybridization solution, including alignment oligo, is made, sample quantified using the nucleic acid feature on the NanoDrop® (Thermo-Fisher) spectrophotometer and the mixture adjusted to ensure optimal concentration for hybridization. Equal amounts of test and reference sample are combined in VWR water; after which the appropriate amount of hybridization solution master mix is added for a total volume of 18 μ L (385K array) or 20 μ L (3x720K array).

The sample/hybridization solution is added to the NimbleGen array/hybridization chamber assembly and incubated with active mixing on the NimbleGen Hybridization System for 40-72 hours at 48°C.

After hybridization is concluded and the array/hybridization chamber assembly is removed from the hybridization system, the hybridization chamber is disassembled and the arrays are carefully washed and dried via centrifugation in a low-light, low-ozone environment followed by immediate scanning.*

Fluorescent Scanning

In the images below, standard resolution (5 μ) scanning and high-resolution (3 μ and 2 μ) scanning was performed using the INNOPSYS InnoScan®900AL per the manufacturer's instructions




Figure 3a. Examples of **1x385K (16 μ spot size)** Roche NimbleGen array images at 5 μ , 3 μ and 2 μ resolution from the INNOPSYS InnoScan®900AL microarray scanner. Note circled areas of corresponding features on each array with successively better resolution as you view the images moving left to right.




Figure 3b. Zoomed-in image of the upper left feature from Circle A in Figure 3a showing increased number of pixels per feature for a 16 μ diameter feature as resolution increases. Note 12 pixels for 5 μ resolution image, 38 pixels for 3 μ image, and 56 pixels for 2 μ images.

Fluorescent Scanning (continued)




Figure 4a. Examples of **3x720K "HD2" (13 μ spot size)** Roche NimbleGen array images at 5 μ , 3 μ and 2 μ resolution from the INNOPSYS InnoScan®900AL microarray scanner. Note circled areas of corresponding features on each array with successively better resolution as you view the images moving left to right




Figure 4b. Zoomed-in image of feature from Circle B in Figure 4a showing increased number of pixels per feature for a 13 μ diameter feature as resolution increases. Note 4 pixels for 5 μ resolution image, 16 pixels for 3 μ image, and 36 pixels for 2 μ images.